

Coriolis Air sampler for bio-contamination






Coriolis[®] μ



www.coriolis-airsampler.com

► Air sampler designed for clean rooms, hospital, indoor and outdoor air control


- High air flow rate: 100 – 300 L/min
- Light: 3 kg
- Collect viable, non cultivable and total flora, pollens, viruses...
- Battery operating (1 hour of sampling)
- Validated according to ISO 14698-1





Air sampling technology

- ▶ **Patented cyclonic technology** concentrating particles from 0.5 to 20 μ m into a sterile liquid collection media
- ▶ Captures and concentrates all **airborne particles** (pollens, spores, bacteria, fungi, virus, allergens, endotoxines...)
- ▶ Qualified by independent agencies : HPA (UK)





Focus on alternative method: cyclonic method

High air flow rate
300 l/min

A liquid sample output




No saturation of the collection media for charged environment and for long time monitoring

A divisible sample for several analysis in parallel
Comparison with the impaction possible

Several types of analyses (PCR, cytometry, immunoanalysis...)
• Shorter response time
• 3 hours protocol validated

Contamination data beyond cultivable flora:
VNC, Pollens, Allergens, Virus...





Rapid detection with Coriolis + RMM

Traditional method (culture): based on cultivability – several days for results

- ▶ Cytometry (ex: Scan RDI) : based on viability – 3 hours
- ▶ ATP Biolum (ex: Pallcheck) : based on viability – 3 to 24 hours
- ▶ ELISA/immunoassay : based on Ag/Ab – max 10 hours
- ▶ PCR : based on DNA or RNA – max 10 hours



- ▶ More reactivity
- ▶ More information
- ▶ More flexibility
- ▶ More expensive in terms of equipment
- ▶ Less expensive in terms of impact on production process



Coriolis consumables

1. Cones with caps (by 50) – 05237-1-003
2. Sterile cones with caps (by 50) – 05237-1-101
3. Collection liquid 15ml doses (by 50) – 05237-1-202



Expiration date : 5 years



Expiration date : 18 months



Expiration date : 9 months




Case Study : Airborne microbiological monitoring in a pharmaceutical production site




Coriolis® - Bertin Technologies

+



Chemscan RDI (Solid phase cytometry) – AES Chemunex


VS.



Impaction Method



22 cleanrooms in a pharma production site (B, C and D grades)



Coriolis® μ + ScanRDI® allowed a better representativeness of the airborne contamination. For D grade room which are the most contaminated rooms, Coriolis® provided an access to more exploitable results.



Mastering innovation



Long time monitoring option for Coriolis μ

Bertin Technologies, a company of **ENIM** Group



NEW in 2013

Long Time Monitoring option for Coriolis μ



- ▶ Bioaerosol samplings with Coriolis μ are limited to 10 minutes as evaporation naturally occurs during collection
- ▶ The long time monitoring option maintains a constant volume of collection liquid over the entire sampling period
- ▶ Main features
 - ▶ From 1 minute to 6 hours of collection
 - ▶ Automatic liquid injection (manually set)
 - ▶ Compatible with any Coriolis μ
- ▶ Provide new possibilities of investigation



Mastering innovation




Compressed Gas Analysis with Coriolis Technology

Bertin Technologies, a company of **ENIM** Group



Coriolis compressed gas description prototype




View of the back



How does this work

► Automation of the collection

- From 10 to 10,000 liters of compressed nitrogen
- For 1 cubic meter, the process is approx 15 minutes







Wet cyclone technology

▶ Coriolis system for particle collection

- ▶ 40 to 70 l/min
- ▶ 10mL of liquid automatically refilled during process (level laser detection)





Specifications of prototype

Use parameters	
Flow rate of gas/air entrance into the Coriolis cone	40 to 70 l/min
Max pressure of gas/air entrance into the Coriolis cone	1 Bar
Sampling volume	10 – 10000 liters
Technical specifications	
Alimentation	100 - 230 V, 0,1A, 50 - 60 Hz
Consumption	20VA
Size & weight	
Large	300mm
Depth	250mm
Total high	350mm
Weight	5Kg
Use conditions	
Temperature	5-40°C
Gas temperature	5-40°C
Humidity	10-90 %
User interface	
Buttons	4 buttons: Valid / + / - / Cancel
Screen	LCD screen, 4 lines



Coriolis: the ideal tool for quality control


SAMPLING



- ✓ Short time of collection (<5 min)
- ✓ Liquid sampling
- ✓ Validated efficiency

+

RMM



- ✓ Culture techniques
- ✓ RMM
- ✓ PCR real time
- ✓ LAL assay
- ✓ Cytometry
- ✓ ADNr16s fingerprint
- ✓ Microarray
- ✓ Microscopic obs.

=

SOLUTION

- ✓ Rapid results
= rapid action
- ✓ Complete information
= more chance of detection

Mastering contamination control