

Western Blot Processors, Shakers, Rockers

WIS BIOMED

1-866-692-1249
www.wisbiomed.com



BlotCycler 6 tray Western Blot Processor

BlotCycler™ automates the process of Western blotting. All steps after protein transfer to a membrane are fully automated: shaking, blocking, hybridization, washing, and recycling of the primary antibody.

- No vacuum, pressurized air or gas, or additional pumps required.
- Can recover primary antibodies.
- Easy programming from the front panel. Program up to 100 protocol steps.
- 6 trays. 6-12 blots per run.
- Audible alarm at the end of the protocol.
- Cleaning cycle at the end of the protocol.



Automated West Blot Processor, Freedom Rocker

Freedom Rocker is an flexible programmable rocker with built-in pumps and reagents delivery and drain machnism. Simply program the rocker and walk away. Freedom Rocker will incubate shake, change buffer and compete your western blot process.

- Freedom Rocker DX4 MD automated western blot gel processor washer
- Programmable Rocker with automatic buffer exchanger, reagent injector
- Designed to process four gels or filters at the time.
- Programmable for rocker speed, time
- Programmable for buffer exchange, recovery
- Fully adjustable reagent injection, recovery



OP-IS50: Incubating mini Orbital Shakers

- OP-IM01: Incubating Microplate Shaker
- OP-OM01: Microplate Shaker
- OP-OS01: Orbital Shaker

- BS-BT4000: Orbi shaker CO2 all environment orbital shaker
- BS-BR3000: Orbi-Shaker Orbital Shaker
- BS-BR5000S: Everlast heavy duty variable speed rocker with stacker
- BS-BR5000: Everlast Heavy duty variable speed rocker



BS-B4000: BioClave 16 liter Autoclave Sterilizer

- Fully automatic; just press start
- Rapid, complete sterilization
- Extremely compact
- Built in water tank, does not require an external water supply

- BS-B1000: BactiZapper MicroSterilizer Cinerator
- an asbestors-free core element that uses infrared heating to produce 815C (1500F) for sterilizing bacterial culture devices



WISBIOMED LLC
120 Boothbay Ave
San mateo, CA 94404

Tel: 1-650-242-8368
Fax: 1-650-242-8244
web: info@wisbiomed.com