



BenchMixer XL

multi-tube vortexer

Hands-free, high capacity vortexing

Up to 2500rpm, a true vortex

Accessories for a variety of tubes,
flasks & plates

Pulsing program for interval mixing

Speed Range: 500 to 2500rpm

Capacity: 2 x 250ml

Orbit: 3.6mm

Timer: 1min. to 99 hr 59min.

Platform Dimensions: 12 x 7 in. / 30 x 17.8 cm

Dimensions: 15 x 9.2 x 16.5 in.

Weight: 13kg / 29lb



50 x 1.5/2.0ml



15 x 50ml



50 x 5, 7 or 10ml



3 x Microplates/96 Well Tube Racks



6 x 50, 100 or 250ml
Volumetric Flasks



9 x 50ml Horizontal
QuEChERS Method