


## Automated Western Blot Technology a better way to do western blot

### how it works:


program the system

add your antibodies,  
blocking and  
wash buffers

walk away



Use the touch screen interface to  
select and design your protocol



Fill tank with washing buffer and  
add up to six different antibodies

**Blot Cycler for four or six trays**  
**Primary antibody Recoverable**  
**Cold room safe operation**

**New touch screen programming**  
**Priced <\$7K, 1yr warranty**  
**Demonstration available**

**BlotCycler ensures high quality and reproducible results for your research**



### FULLY AUTOMATED

- process multiple blots
- all blotting solutions delivered automatically
- hybridize six different antibodies at the same time
- easy programming of all hybridization steps for multiple users



### RELIABLE & EFFICIENT

- reliable detection of low expressing proteins
- highly reproducible results
- precise delivery of all reagents
- automatic primary antibody collection



### USER FRIENDLY

- automatic cleaning
- minimal bench space required
- guaranteed to save time & money
- energy efficient

visit [wisbiomed.com](http://wisbiomed.com) for more information