

Polyethylene Vial Sets

For high throughput processing of seeds, leaves and animal tissues

WISBIOMED

120 Boothbay Ave.
San Mateo CA, 94404



TEL (650) 242-8368

FAX (650) 242-8244

Email: info@wisbiomed.com

Sample Preparation Products

24 Well PE Vial Set (Optional Lined Caps)



Polyethylene vs. Polycarbonate vial

For labs using solvents incompatible with polycarbonate, the 24 Well PE Vial Sets are a useful alternative. Resistant to many solvents that melt polycarbonate, these cost effective PE Vial Sets are suitable for most small seeds and animal tissue. These are not optically clear and not recommended for corn kernels.

Product Number	Description	Quantity	Cost
OP-PEVS 04-240-04	24 Well PE Vial Sets	Case of 10	\$202
OP-PEVS 04-240-05	24 Well PE Vial Sets with lined Caps	Case of 10	\$223

Related Products

OP-PCVS 15-050-01	15 ml PC Vial Sets	Case of 10	\$249
OP-PCVS 04-240-02	24 Well PC Vial Sets	Case of 10	\$249
OP-PCVS 04-240-03	24 Well PC Vial Sets with lined Caps	Case of 10	\$270

WISBIOMED

120 Boothbay Ave.
San Mateo CA 94404