

Design Comparison: Precellys 24 vs. FastPrep 24



Metal Sample Holder

**Metal Sample holder
vs.
Plastic holder**




Automatic Pressured Tube Securing Mechanism

**Vacuum Pressured
vs.
Screw cap**



Completely sealed sample chamber, easy to clean

**Containment seal
vs.
Open air gap**



More efficient engine cooling

**Vent Fan & Grid
Large opening & Metal
vs.
Small, plastic**



Clean Efficient Sample Cooling

**Insulated Lid, Liquid
nitrogen or dry ice chamber
vs.
entire shake unit**



visit wisbiomed.com for more information