

Precellys Emulsion Kit

Get Perfect Emulsions



- **Easy:** add your Adjuvant and Antigen solutions in the tube and let the instrument mix them for you
- **Fast & Efficient:** reproducible high-quality emulsions in only 2 min on Minilys or Precellys Evolution
- **Flexible:** prepare 2 to 7mL of emulsion
- **Cost-effective:** no waste of expensive reagents
- **Safe:** no risk of accidental exposure or Repetitive Stress Injuries (RSI)

An efficient method...

With the same volume of reagents (7 mL)


In 2 minutes




In 15-20 minutes

**Precellys Emulsion kit
with Minilys**

Two Syringe Method



70 injections



Confidence
in reproducibility

50 injections






Risk of Repetitive
Stress Injuries
(RSI)

...to prepare homogenous and stable emulsions

Precellys
Emulsion kit

2-syringe method

Vortex



VISIT OUR WEBSITE

www.bertin-instruments.com

Bertin Technologies - Head Office

10 bis Avenue Ampère - Parc d'activités du Pas du Lac
78180 Montigny-le-Bretonneux FRANCE