



# ***X-CLARITY™ Tissue Clearing System***

***Easy installation, ultrafast clearing, and reliable results***



Cleared Whole Mouse Brain  
ETC 24 hours with X-Clarity™



Cleared Mouse Spinal Cord  
ETC 2 hours with X-Clarity™

The **X-CLARITY™** is the first commercial turn-key system that is optimized to perform the groundbreaking CLARITY method and will greatly simplify the CLARITY tissue clearing for your lab. Many researchers have been unsuccessful in replicating the CLARITY method due to its complex, do-it-yourself procedure. The X-CLARITY ETC eliminates the need to build your own device and has a simple to perform work flow and high rate of success for end users. X-CLARITY incorporated many improvements to the original ETC process.

## **Key Features:**

- Timed polarity switching electrophoresis (Patent pending)
- Fully integrated temperature control with built-in temperature probe and buffer re-circulating pump
- ETC chamber with optimized buffer flow and plate electrodes

## **Want more detailed information about X-CLARITY™ Tissue Clearing System**

For details, visit our website at [www.logosbio.com](http://www.logosbio.com) or contact us at [info@logosbio.com](mailto:info@logosbio.com) or [kim65@logosbio.com](mailto:kim65@logosbio.com).

For in lab demonstration, contacting us at (703) 622-4660 or (703) 942-8867.