

Rocker with automated buffer exchanger



It's a rocker with a pump and a timer!

*The Freedom Rocker™ is a **very quiet** laboratory platform rocker with a pump and a timer. It automatically processes blots and gels. It frees you from babysitting your samples by relieving you from the frequent pouring of liquids into and out of trays.*

- increase productivity
- reduce downtime
- improve safety
- 30 days full refund policy, risk free

> Increase Productivity

The Freedom Rocker™ automates the processing of blots and gels, freeing up time for more demanding activities. Because it automatically dispenses and aspirates reagent from multiple trays simultaneously, the samples in all the trays will experience the exact same protocol. Through increased productivity, the Freedom Rocker™ pays for itself in about 30 uses.

> Reduce Downtime and Shorten Cycles

The Freedom Rocker™ allows you to process blots and gels anytime . . . even during meetings or overnight. No need to wait until the next morning to start processing.

> Improve Safety

The Freedom Rocker™ automatically drains hazardous reagents from trays into the proper containers, reducing exposure to dangerous materials.

> Risk Free Purchase

The Freedom Rocker™ comes with a 30 day money back guarantee and a 2 year warranty.

Order Instruction:

Fax or email Purchase Order to WISBiomed at 650-242-8244, sales@wisbiomed.com

Question: call 650-242-8368

BASIC VERSION FEATURES

- Compact. Uses less than 10 in. of lab bench. 9.8 in. (25 cm) wide by 13.7 in. (35 cm) deep, including the rocker. (12 inches wide for the 4 tray version)
- No vacuum, pressurized air or gas, or additional pumps required.
- Rocker platform 9.8 in. (25 cm) wide (FRSK22) or 12 in. (30 cm) wide (FRSK42) by 12 inch deep (30 cm) rocker platform standard 16 in. (41 cm) wide by 15 in. (38 cm) deep platform available.
- Can recover reagent solution.
- Unattended dispensing and aspiration of liquids to 1 or 2 trays (FRSK22 or FRSK42) or 4 trays (FRSK42).
- Four adjustable Parameters: rocking speed, reagent time, wash buffer time, dispense volume.
- Rocking continues after the final rinse and the wash buffer remains in the trays to keep the blots or gels wet.
- Warranty: 2 years parts and labor.
- 30 day money back guarantee.

PROGRAMMABLE (DX) VERSION

- All of the above features
- Can program up to 9 protocol steps
- Clearing rinse can be activated
- Alarm at the end of the last step can be activated or deactivated.
- Color coded tubing keep connections to bottles simple.