


Eyepiece Digital Camera




Eyepiece Digital Camera

- Eyepiece Adapter Digital Camera
- USB 2.0 Output
- 1.3M / Resolution 1280x 1024
- 70mm x 45mm (D)
- 36mm inside diameter

Work with most microscopes
Just attach to you eyepiece
Save time, Capture Images
Easy to share and teach

1.3MP Digital Microscope-USB ready



- Handheld Microscope
- 10X~50x,200x Continuous Mag's
- digital "zoom in" VGA up to 330x
- 1.3M / Resolution 1280x 1024
- * AM413T w/Measurement feature
(not available on AM411T or Mac OS)
- Built-in 8 white LED's

Many models are available: handheld, white light, UV-White light, long view, polarize kit