

# Benchmark Hotplates, Stirrers & Hotplate Stirrers



The Benchmark hotplates, stirrers and hotplate-stirrers feature an exceptionally durable, chemical resistant white ceramic work surface. Their space-efficient design (8 X 9 in. footprint), makes them ideal for use on crowded benchtops and inside of bio-hoods.

Advanced microprocessor controls with convenient turn knobs allow quick, precise adjustment and maintenance of speed and temperature. Safety indicator LED's on the front panel indicate when the heating and/or stirring functions have been activated.

With a square 7.5 in work surface, all three models are compatible with a wide variety of popular sizes of borosilicate glass beakers, flasks, bottles and other vessels. The hotplate/magnetic stirrer model (H4000-HS) comes complete with a support rod for mounting thermometers and temperature probes.

- + Large, white ceramic top plate
- + Compact, 8 x 9" footprint
- + Three models:  
Hotplate, Stirrer and Hotplate Stirrer

## Technical Data

### **Hotplate (H4000-H)**

Temperature Range: Amb. +5 to 380 °C\*  
Platform: 7.5 x 7.5 in./ 19 x 19 cm  
Control: Quick adjustment knob  
Dimensions: 8 (W) x 9 (D) x 4.5 (H) in.  
Electrical Data: 120V, 60 Hz / 230V, 50/60 Hz

### **Magnetic Stirrer (H4000-S)**

Speed Range: 60 to 1500 rpm  
Platform: 7.5 x 7.5 in./ 19 x 19 cm  
Control: Quick adjustment knob  
Dimensions: 8 (W) x 9 (D) x 4.5 (H) in.  
Electrical Data: 120V, 60 Hz / 230V, 50/60 Hz

### **Hotplate Stirrer (H4000-HS)**

Temperature Range: Amb. +5 to 380 °C\*  
Speed Range: 60 to 1500 rpm  
Platform: 7.5 x 7.5 in./ 19 x 19 cm  
Control: Quick adjustment knob  
Dimensions: 8 (W) x 9 (D) x 4.5 (H) in.  
Electrical Data: 120V, 60 Hz / 230V, 50/60 Hz

## Ordering Information

H4000-H Benchmark hotplate, 115V\*  
H4000-S Benchmark magnetic stirrer, 115V\*  
H4000-HS Benchmark hotplate magnetic stirrer, 115V\*

\* To order a product in 230V, please add -E to item number.

**Distributed by WISBIOMED LLC**  
**1-866-692-1249 wisbiomed.com**