

myFUGETM Mini

Mini Centrifuge



**Eight Tube
Capacity**

Benchmark Scientific

The new MyFUGE Mini centrifuge is a personal centrifuge that virtually fits in the palm of your hand. Ideal for quick spin downs of microtubes and PCR tubes, the MyFuge is extremely easy to use. Simply close the lid and the rotor instantly reaches 6,000 rpm. Open the lid, and the rotor quickly decelerates to a stop for removing samples.

An increased capacity, eight-position microtube rotor is included and comes installed. A second rotor for 0.2ml PCR tubes and strips is also included and stored on the bottom of the centrifuge. When the rotor is not in use, storage in this compartment prevents loss or misplacement of the rotor.

- + Eight position microtube rotor
- + Storage compartment for PCR rotor
- + Near silent operation
- + Starts and stops in just seconds
- + Conserves valuable bench space

Technical Data

Speed:	6,000rpm / 2,000 xg
Capacity:	8 x 1.5 / 2.0ml tubes 16 x 0.2ml PCR tubes 2 x PCR strips (8x0.2ml)
Dimensions:	4.5 x 5.9 x 4.5 in 11.4 x 15 x 11.4 cm
Weight:	1.1 lbs / 5 kg
Electrical:	100 - 240V, 50-60Hz



Storage Compartment



Includes Both Rotors

Ordering Information

MyFUGE Mini Centrifuge, blue lid
MyFUGE Mini Centrifuge, clear lid
MyFUGE Mini Centrifuge, neon green lid
MyFUGE Mini Centrifuge, purple lid
MyFUGE Mini Centrifuge, red lid

*All include microtube rotor, PCR rotor, 0.5ml & 0.2ml adapters

Distributed by WISBIOMED LLC
1-866-692-1249 wisbiomed.com