

myBLOCK™

Mini DryBath



- + Compact, fits in the palm of your hand
- + Exchangeable blocks, for tubes 0.2 to 50ml
- + Plastic cover ensures temp. uniformity
- + Digital temperature control
- + Portable, 12V vehicle power adapter

Technical Data

Temp. Range:	Ambient +5 to 80°C
Temp. Accuracy:	+/- 0.5°C
Temp. Increments:	0.1°C
Temp. Uniformity:	+/- 0.2°C
Timer:	1 to 99 min.
Block Construction:	High Grade Aluminum
Dimensions:	4.4 x 5.9 x 4 in. 11.2 x 15 x 10 cm
Weight:	1.1 lbs / 5 kg
Electrical:	100 - 240V, 50-60Hz*

* 12V vehicle power adapter also available



Ordering Information

BSH100 MyBlock Mini Digital Dry Bath with block lifter (Blocks sold separately)

BSH100-02	Block, 40 x 0.2ml (or X x 0.2ml strips)
BSH100-05	Block, 24 x 0.5ml
BSH100-15	Block, 15 x 1.5ml
BSH100-1520	Block, 15 x 1.5 or 2.0ml
BSH100-150	Block, 4 x 15ml*
BSH100-500	Block, 2 x 50ml*
BSH100-A12	12V vehicle power adapter

*Not compatible with plastic cover

Distributed by WISBIOMED LLC
1-866-692-1249 wisbiomed.com