

Automated Rocker with buffer exchanger



fully programmable rocker

enables you to be productive while
processing gels, membranes, or slides

MANUAL STEPS

POUR SOLUTION/
REAGENT INTO TRAY
AND START ROCKER

POUR OUT
SOLUTION

POUR IN
WASH BUFFER

POUR OUT
WASH BUFFER

POUR IN
WASH BUFFER

POUR OUT
WASH BUFFER

POUR IN
WASH BUFFER

POUR OUT WASH
BUFFER

AUTOMATED STEPS

POUR SOLUTION/
REAGENT INTO TRAY
AND PRESS **'RUN'**

**YOU ARE
FREE TO
WORK
ON OTHER
TASKS**

**FINISH AUTO
CYCLE**

For processing of DNA
gels, protein gels, and
Western blots

Supports 3 buffers
or reagents
Automatic buffer
filling, draining
reagent injection
and recovery

Programmable
rocker with 2 or
4 gel trays

Basic version has
limited programs

Deluxe version has
9 programmable steps

Package includes
2-year warranty

Call 650-242-8368
for quotations or
visit WISBIOMED.com