


myBATH™ 4L

The complete personal water bath



- + Hinged lid is also removable
- + Includes tube racks for 1.5/2.0ml, 15 and 50ml tubes
- + Quik-Cal, in-lab user calibration function
- + Digital accuracy and convenience

Included Tube Racks:



This compact 4 liter water bath features easy to use controls that allow the user to digitally select and monitor temperature. Its hinged lid provides a secure, covered environment and can be flipped open (to the rear) or completely and instantly removed at the user's option.

The MyBath is also the only water bath that comes with test tube racks for 50ml, 15ml and 1.5/2.0ml tubes, making it complete and ready to use with all of the most frequently used tube sizes.

Designed with the molecular biology lab in mind, the MyBath 4 liter is the complete personal water bath.

Technical Data

Temp. Range:	Ambient +5 to 100°C
Capacity:	4L
Temp. Accuracy:	+/- 0.2°C
Temp. Increments:	0.1°C
Temp. Uniformity:	+/- 0.2°C
Exterior Dimensions:	7.7 x 13.9 x 9 in. / Height: 12.2 in. with lid
(WxDxH)	19.5 x 35 x 23 cm / Height: 31 cm with lid
Chamber Dimensions:	5.3 x 9.4 x 5.5 in. / 13.5 x 24 x 14 cm
Weight:	15 lbs / 6.5 kg
Electrical:	115 or 230V, 50-60Hz, 400W

Ordering Information

B2000-4* MyBath 4L Digital Water Bath, includes tube racks for 1.5/2.0ml, 15 & 50ml tubes

B2000-T5 Optional tube rack for 40 x 0.5ml tubes

* 115V with US plug. For 230V (EU plug), please add (-E)

Distributed by WISBIOMED LLC
1-866-692-1249 wisbiomed.com